

# OrthoPulse™ Ultra 50

Compact Advanced Radial Pulse Wave  
System for the Treatment of Acute and  
Chronic Musculoskeletal Pain

Extracorporeal Pulse Activation  
Technology (EPAT®) System


Powered by STORZ Medical

Regenerative Medicine Redefined

STORZ MEDICAL

CuraMedix  
Innovative Technologies for Advanced Healing

OrthoPulse™ Ultra 50





## The OrthoPulse™ Ultra 50

The OrthoPulse™ Ultra 50 acoustic pressure wave device is a light weight, portable EPAT device ideal for medical professionals in private practice or in clinics; available at an affordable price with low maintenance costs.

- Facts & Figures**
- Built-in compressor
  - Weight: 21 lbs.
  - Application pressure: 0.3-4.0 bar
  - Control and display all on handpiece
  - Pulse frequency: 1-17 Hz
  - Vibration therapy (V-ACTOR® Handpiece)

Widely used in podiatry, chiropractic, sports medicine, pain management, physical therapy in both the U.S. and around the world, this evidence-based, non-invasive approach is extremely effective in treating the root cause of the pain and increasing mobility.

### Applications

Generally, acute or chronic musculoskeletal pain and/or pain that significantly impairs mobility or quality of life including:

- Foot and ankle
- Lower extremity
- Upper extremity
- Back and chest
- Shoulder
- Neck
- Myofascial trigger points



### The Design

The pioneering Single Frame Casing and innovative three-dimensional design of the OrthoPulse™ Ultra combines unmatched performance, reliability and value.

### The Handpiece

System functionality and Individual Parameter Settings (IPS-Controls) have been integrated into the new "Active-Tip-Control" display on the handpiece make treating patients easier.

### The CuraMedix Partner Portal

Access to the latest regenerative medicine treatment information including research, articles, physician and patient testimonials, and more.

**Patient Care Redefined**  
  
Innovative Technologies for Advanced Healing

TOLL FREE: 877-699-8399 TEL: 401-333-6500 FAX: 401-633-6565  
 INFO@CURAMEDIX.COM • WWW.CURAMEDIX.COM

#### Important Information for U.S. Customers:

Certain devices and references made herein to specific indications of use may have not received clearance or approval by the United States Food and Drug Administration. Practitioners in the United States should first consult with their local CuraMedix representative in order to ascertain product availability and specific labeling claims. Federal (USA) law restricts certain devices referenced herein to sale, distribution, and use by, or on the order of a physician, dentist, veterinarian, or other practitioner licensed by the law of the State in which she/he practices to use or order the use of the device.